

Multiwell Plates

(6 - 12 - 24 - 48 well)

For both adherent and suspension cell culture, choose Greiner Bio-One's Multiwell plates in either 6, 12, 24 or 48 well format. The special chimney wells in conjunction with low evaporation lids significantly reduces the risk of cross contamination between the wells. The plates are easily stacked, both with and without lids. Due to their advanced design, the plates still exhibit important features like constant heat- and gas exchange even when stacked. Uniform plate dimensions (well depth, -diameter and distance between the wells) comply with SBS standards, enabling easy manual and automated handling. Each plate is sterile, non-pyrogenic, DNase/RNase-free, and has a printed lot number for traceability.



Features

Plate format	6 / 12 / 24 / 48 wells
Growth area	9.6 / 3.9 / 1.9 / 1.0 cm ²
Working vol. [ml]	2.5 / 2.4 / 0.5-1.5 / 0.5-1
<ul style="list-style-type: none"> ■ Low evaporation lid ■ Chimney well design ■ Alpha-numeric well coding ■ Writing area and lot no. on each plate ■ Easily stackable ■ Non-pyrogenic, RNase- & DNase-free 	

Ordering Information

Cat.-No.	Description	Qty/bag	Cat.-No.
657 160	6 well plate, tissue culture treated, with lid ,sterile	1	657 160
665 180	12 well plate, tissue culture treated, with lid, sterile	1	665 180
662 160	24 well plate, tissue culture treated, with lid, sterile	1	662 160
677 180	48 well plate, tissue culture treated, with lid, sterile	1	677 180
657 185	6 well plate, suspension culture quality, with lid, sterile	1	657 185
665 102	12 well plate, suspension culture quality, with lid, sterile	1	665102
662 102	24 well plate, suspension culture quality, with lid, sterile	1	662 102
677 102	48 well plate, suspension culture quality, with lid, sterile	1	677 102

Revision: August 2003 – 073 025

www.gbo.com/bioscience

Deutschland (Zentrale)
Greiner Bio-One GmbH
Tel: (+49) 7022 948-0
E-Mail: info@de.gbo.com

Belgien
Greiner Bio-One N. V.
Tel: (+32) 2-4 61 09 10
E-Mail: info@be.gbo.com

Frankreich
Greiner Bio-One SAS
Tel: (+33) 169-86 25 50
E-Mail: infos@fr.gbo.com

Japan
Greiner Bio-One Co. Ltd.
Tel: (+81) 3-35 05-88 75
E-Mail: info@jp.gbo.com

Österreich
Greiner Bio-One GmbH
Tel: (+43) 7583 6791-0
E-Mail: office@at.gbo.com

Niederlande
Greiner Bio-One B. V.
Tel: (+31) 172-42 09 00
E-Mail: info@nl.gbo.com

UK
Greiner Bio-One Ltd.
Tel: (+44) 1453-82 52 55
E-Mail: info@uk.gbo.com

USA
Greiner Bio-One North America Inc.
Tel: (+1) 704-261-78 00
E-Mail: info@us.gbo.com


greiner bio-one