

CELLMASTER™

Roller Bottles with Filter Cap

Greiner Bio-One Polystyrene Roller Bottles with a smooth or ribbed surface range in size from 850 cm² to 4720 cm² growth area (resp. 270 or 500 mm height) and are available either with standard or filter screw caps. The filter membrane is a hydrophobic capillary pore membrane with a defined pore size of 0.2 µm to guarantee optimal protection against contamination while providing a high gaseous exchange rate which enables a mass cultivation of cells using bicarbonate buffer and CO₂.



Features

- High grade polystyrene
- "Best Before" date and lot number on each roller bottle
- Sterile
- Non-pyrogenic (limit <0,06EU/ml)
- DNase-/RNase- free
- Tested by USP Class VI
- Filter Cap with capillary pore membrane
- Defined pore size 0.2 µm
- Excellent gas exchange

Ordering Information

Cat.-No.	Roller Bottle with Filter Cap	Qty/bag	Qty/case
680 045	1 X PS roller bottle, 850 cm ² , smooth surface, filter cap, 122 x 270 mm	24	48
680 058	1 X PS roller bottle, 850 cm ² , smooth surface, filter cap, 122 x 270 mm	2	24
680 068	1 X PS roller bottle, 850 cm ² , smooth surface, filter cap, 122 x 270 mm	24	24
681 062	2 X PS roller bottle, 1700 cm ² , expanded surface, filter cap, 122 x 270 mm	2	24
681 072	2.5 X PS roller bottle, 2125 cm ² , expanded surface, filter cap, 122 x 270 mm	2	24
682 015	1 XL PS roller bottle, 1700 cm ² , smooth surface, filter cap, 122 x 500 mm	12	12
682 025	1 XL PS roller bottle, 1700 cm ² , smooth surface, filter cap, 122 x 500 mm	12	24
682 065	1 XL PS roller bottle, 1700 cm ² , smooth surface, filter cap, 122 x 500 mm	1	12
682 078	5 XL PS roller bottle, 4720 cm ² , expanded surface, filter cap, 122 x 500 mm	1	12
383 382	Filter Screw Cap, single packed, sterile	1	50

Revision: December, 2006

www.gbo.com/bioscience

Germany (Main office)
Greiner Bio-One GmbH
Phone: (+49) 7022 948-0
E-Mail: info@de.gbo.com

Belgium
Greiner Bio-One N. V.
Phone: (+32) 2-4 61 09 10
E-Mail: info@be.gbo.com

France
Greiner Bio-One SAS
Phone: (+33) 1 69-86 25 50
E-Mail: infos@fr.gbo.com

Japan
Greiner Bio-One Co. Ltd.
Phone: (+81) 3-35 05-88 75
E-Mail: info@jp.gbo.com

Austria
Greiner Bio-One GmbH
Phone: (+43) 7583 6791-0
E-Mail: office@at.gbo.com

Netherlands
Greiner Bio-One B. V.
Phone: (+31) 172-42 09 00
E-Mail: info@nl.gbo.com

UK
Greiner Bio-One Ltd.
Phone: (+44) 1453-82 52 55
E-Mail: info@uk.gbo.com

USA
Greiner Bio-One North America Inc.
Phone: (+1) 704-261-78 00
E-Mail: info@us.gbo.com



greiner bio-one