

# Media Bottle

Tetragonal media bottle made of ultra-clear polyethylene terephthalate (PET).

Media bottle with screw cap for liquid storage. Thermal stability from -60°C up to +60°C.

- Available in three sizes: 100, 500 and 1000ml
- With graduation
- Sterile
- Free of DNase/RNase
- Non-cytotoxic and non-pyrogenic



## Height x Width in mm

- 100ml-flask: 108.3 x 53.0
- 500ml-flask: 177.3 x 77.0
- 1000ml-flask: 218.8 x 92.0

## Ordering Information

Cat.-No.	Media Bottle	Qty/tray	Qty/case
951 700	Media Bottle 100ml	25	100
950 700	Media Bottle 500ml	25	50
952 700	Media Bottle 1000ml	12	24

Revision: Juli 2006

[www.gbo.com/bioscience](http://www.gbo.com/bioscience)

**Germany (Main office)**  
Greiner Bio-One GmbH  
Phone: (+49) 7022 948-0  
E-Mail: [info@de.gbo.com](mailto:info@de.gbo.com)

**Belgium**  
Greiner Bio-One N. V.  
Phone: (+32) 2-4 61 09 10  
E-Mail: [info@be.gbo.com](mailto:info@be.gbo.com)

**France**  
Greiner Bio-One SAS  
Phone: (+33) 169-86 25 50  
E-Mail: [infos@fr.gbo.com](mailto:infos@fr.gbo.com)

**Japan**  
Greiner Bio-One Co. Ltd.  
Phone: (+81) 3-35 05-88 75  
E-Mail: [info@jp.gbo.com](mailto:info@jp.gbo.com)

**Austria**  
Greiner Bio-One GmbH  
Phone: (+43) 7583 6791-0  
E-Mail: [office@at.gbo.com](mailto:office@at.gbo.com)

**Netherlands**  
Greiner Bio-One B. V.  
Phone: (+31) 172-42 09 00  
E-Mail: [info@nl.gbo.com](mailto:info@nl.gbo.com)

**UK**  
Greiner Bio-One Ltd.  
Phone: (+44) 1453-82 52 55  
E-Mail: [info@uk.gbo.com](mailto:info@uk.gbo.com)

**USA**  
Greiner Bio-One North America Inc.  
Phone: (+1) 800-884-47 03  
E-Mail: [info@us.gbo.com](mailto:info@us.gbo.com)

  
**greiner bio-one**